


Peripheral IV connector range:



Chalfont Connector



Coventry Connector



High Flow Y Connector



2-way Safe-Connect



Ordering information

Product Code	Product Description	Box Quantity
6507	Chalfont Connector	50
6510	Coventry Connector	50
6517	High Flow Y Connector	50
6700	2-way Safe-Connect	50

Clinical data

¹Safe Anaesthesia Liaison Group (SALG), Guaranteeing Drug Delivery in Total Intravenous Anaesthesia, October 2009

Mediplus are experts in the delivery of high performance medical devices that support positive patient outcomes in the following fields -



Mediplus Ltd.
 +44 (0) 1494 551200
 help@mediplus.co.uk
 The Gateway Centre, Coronation Road,
 High Wycombe HP12 3SU
 www.mediplus.co.uk



INVESTORS IN PEOPLE™
 We invest in people Silver



Mediplus Peripheral IV Connectors

Safer, more efficient drug delivery

Mediplus Peripheral IV Connectors

Safer, more efficient drug delivery

Mediplus peripheral IV connectors are designed to ensure optimal drug delivery for enhanced safety in ICU, HDU, orthopaedic and obstetric departments. These high-performance connectors offer a wide range of benefits including:

Enables rapid emergency response – the high-flow anti-reflux valve (Coventry Connector and High Flow Y Connector) offers a uniquely fast flow rate for rapid IV fluid delivery and prevents back flow into the IV line

Greater control – the integrated anti-siphon valve prevents free flow and back flow, reducing the risk of over or under dosing during infusion

Comfort & convenience – the removable bridge (Coventry, 2-way Safe-Connect and Chalfont Connectors) helps prevent the device becoming dislodged and keeps it compact for greater patient comfort and manoeuvrability

Ergonomic design – the J-loop design (Coventry, 2-way Safe-Connect and Chalfont Connectors) prevents kinking and occlusion of tubing for unimpeded flow

Reduced risk – the needle-free access ports (2-way Safe-Connect) support infection control and reduce the risk of needle stick injuries

Wide range – single, double or triple infusion options available

“I have been using the Coventry Connectors for more than 10 years now and I find them really useful. The multiple ports make it easy to run multiple infusions and fluids together. They are kink resistant and the rate of fluid flow is really good”

Payal Kajekar,
Luton & Dunstable Hospital

Key features



J-loop to prevent kinking/occlusion of the tubing



Removable bridge offers more flexibility when moving patients



Anti-siphon valves to prevent back flow and free flow



Wide bore gravity line to allow high flow of IV fluid



Anti-reflux valve to prevent back flow



Needle-free access port

